



## Get better with **ActivaScrew™ Interference**

High Usability for Various Indications



Choose


**ActivaScrew™**  
**Interference**

for its unique material properties.

Hotline 0800 588 8005  
[www.bluerockmedical.de](http://www.bluerockmedical.de)

# Get better with ActivaScrew™ Interference

High Usability for Various Indications



## 1. Hand

- Thumb Caprometacarpal (CMC) Ligament Reconstruction

## 2. Elbow

- Medial Collateral Ligament (MCL) Reconstruction
- Lateral Collateral Ligament (LCL) Reconstruction

## 3. Shoulder

- Biceps Tenodesis

## 4. Knee

- Anterior Cruciate Ligament (ACL) Reconstruction
- Posterior Cruciate Ligament (PCL) Reconstruction
- Medial Patellofemoral Ligament (MPFL) Reconstruction

## 5. Foot and ankle

- Lateral Ankle Reconstruction
- Flexor Digitorum Longus (FDL) Tendon Transfer
- Flexor Hallucis Longus (FHL) Tendon Transfer
- Achilles Tendon
- ActivaScrew

The ActivaScrew™ Interference bioabsorbable screws are constructed of bioabsorbable poly(L-lactide-co-glycolide) (PLGA).

ActivaScrew™ Interference

Patents: FI 124190. EP, US, FI Patent Pending

# ActivaScrew™ Interference Offer

## HIGH USABILITY FOR VARIOUS INDICATIONS

ActivaScrew™ Interference is indicated for fixation of tissue, including ligament or tendon to bone, or a bone-tendon to bone in the presence of appropriate immobilization or controlled mobilization. Interference fixation is appropriate for surgeries of the knee, shoulder, elbow, ankle, foot and hand/wrist.

## WIDE SELECTION OF DIFFERENT SIZES

Diameters of 4 - 10 mm and lengths from 10 - 33 mm.

## EASY INSERTION ALONG THE GUIDE WIRE (K-WIRE)

The ActivaScrew™ Interference is cannulated. It allows insertion of the screw over a guide wire (1.6 or 0.8 mm) facilitating screw placement.

## EASY POSITION IMAGING DURING THE OPERATION

## CENTER DRIVEN TORX DRIVER FOR OPTIMAL INSERTION

## HIGH STRENGTH PROPERTIES OFFER EASY INSERTION AND SAFE MEDICAL USE

High initial mechanical properties and isoelasticity of ActivaScrew™ Interference are created during the manufacturing of the product. These features increase the reliability of the screw in clinical use.

## SELF-LOCKING SL™

After operation, dimensions of ActivaScrew™ Interference slightly change improving screw's fit inside the tunnel and thus reducing the risk of an unstable fixation. The Self-Locking SL™ property is activated under hydrolytic conditions.

## ISOELASTICITY

The bending modulus of the ActivaScrew™ Interference is closer to the value of cortical bone compared to metallic implants.

## MR SAFE

Metallic implants are known to disturb MR imaging. Bioabsorbable materials are not magnetic and do not disturb MR imaging.

## ENABLES DRILLING THROUGH THE PRODUCT IN NEED OF RE-OPERATION

## NO NEED FOR REMOVAL OPERATION

## IMPLANT IS SUPPLIED GAMMA STERILIZED

It is safe, free of gas remnants and the risk of cross infection is reduced.



# Example indications for ActivaScrew™ Interference

## Anterior Cruciate Ligament (ACL) Reconstruction

Bioretec offers a wide selection of different screw sizes and easy to use instrumentation for knee ligament reconstruction. Self-Locking SL™ property and graft friendly design of ActivaScrew™ Interference increase the reliability of the screw in clinical use.

In anterior cruciate ligament (ACL) reconstruction, a cannulated Notcher (B-II-8010) can be used to prepare the insertion point for the screw into the femoral and tibial tunnels. Using a cannulated torx driver for a 7/8/9/10 mm interference screw (B-II-9031) the ActivaScrew™ Interference is picked up from the holder and inserted along the guide wire (B-II-1624B) into the tunnel.

Typically, 7 - 10 mm screws are used for ACL reconstruction.



## Medial Patellofemoral Ligament (MPFL) Reconstruction

In addition to 7 - 10 mm interference screws, the Bioretec product range includes 4 - 6 mm interference screws for knee, shoulder, elbow, ankle, foot, hand and wrist surgery. Screws are cannulated and easy to use with a center driven Torx driver. Isoelasticity, Self-Locking SL™ and high strength properties enable easy insertion, safe medical use and good treatment results.

In case of medial patellofemoral ligament (MPFL) reconstruction a 4 - 5 mm ActivaScrew™ Interference can be used to fix a tendon graft onto the patella and 5 - 6 mm screws to fix a tendon graft onto the femur. 4 - 6 mm interference screws are inserted using a Torx shaft for 4/5/6 mm interference screw (B-II-9032). The Torx shaft is easy to attach to the large handle (B-ISC-9013) with quick coupling.






## Example indications for ActivaScrew™ Interference



### Lateral Ankle Reconstruction

For lateral ankle reconstruction, 6 mm ActivaScrew™ Interference can be used to fix a tendon graft to achieve an anatomic reconstruction of the anterior talofibular and calcaneofibular ligaments in the ankle. A wide selection of different screw diameters and lengths enable you to choose the most suitable screw for the situation.

The ActivaScrew™ Interference delivery concept supports an aseptic way of working during the surgery. The ActivaScrew™ Interference is easy to pick up from the holder using a center driven Torx driver (B-II-9031, B-II-9032). The Torx driver extends through the length of the screw transmitting force throughout the screw making insertion of the screw easy and secure.



The ActivaScrew™ Interference bioabsorbable screws are constructed of bioabsorbable poly(L-lactide-co-glycolide) (PLGA).

**ActivaScrew™ Interference**

Patents: FI 124190. EP, US, FI Patent Pending

# ActivaScrew™ Interference product range



## REMARK

Prior to using Bioretec Implants examine thoroughly the Instructions For Use - inside each product package.

## ORDERING INFORMATION

For more information please contact your nearest distributor or our customer service at [info@bluerockmedical.de](mailto:info@bluerockmedical.de). Valid product range available from our website: [www.bluerockmedical.de](http://www.bluerockmedical.de).

## Product range

Product Code	Description	Product Code	Description
B-AI-0410	ActivaScrew™ Interference 4 x 10 mm	B-AI-0820	ActivaScrew™ Interference 8 x 20 mm
B-AI-0415	ActivaScrew™ Interference 4 x 15 mm	B-AI-0824	ActivaScrew™ Interference 8 x 24 mm
B-AI-0420	ActivaScrew™ Interference 4 x 20 mm	B-AI-0828	ActivaScrew™ Interference 8 x 28 mm
B-AI-0510	ActivaScrew™ Interference 5 x 10 mm	B-AI-0830	ActivaScrew™ Interference 8 x 30 mm
B-AI-0515	ActivaScrew™ Interference 5 x 15 mm	B-AI-0833	ActivaScrew™ Interference 8 x 33 mm
B-AI-0520	ActivaScrew™ Interference 5 x 20 mm	B-AI-0920	ActivaScrew™ Interference 9 x 20 mm
B-AI-0615	ActivaScrew™ Interference 6 x 15 mm	B-AI-0924	ActivaScrew™ Interference 9 x 24 mm
B-AI-0620	ActivaScrew™ Interference 6 x 20 mm	B-AI-0928	ActivaScrew™ Interference 9 x 28 mm
B-AI-0626	ActivaScrew™ Interference 6 x 26 mm	B-AI-0930	ActivaScrew™ Interference 9 x 30 mm
B-AI-0720	ActivaScrew™ Interference 7 x 20 mm	B-AI-0933	ActivaScrew™ Interference 9 x 33 mm
B-AI-0724	ActivaScrew™ Interference 7 x 24 mm	B-AI-1024	ActivaScrew™ Interference 10 x 24 mm
B-AI-0726	ActivaScrew™ Interference 7 x 26 mm	B-AI-1028	ActivaScrew™ Interference 10 x 28 mm
B-AI-0730	ActivaScrew™ Interference 7 x 30 mm	B-AI-1030	ActivaScrew™ Interference 10 x 30 mm

## ActivaScrew™ Interference Instrumentsa

Product Code	Description
B-II-9032	Torx shaft for 4/5/6 mm interference screw with quick coupling
B-ISC-9013	Large handle, cannulated, quick coupling
B-II-9031	Torx driver for 7/8/9/10 mm interference screw
B-II-8010	Notcher for 7/8/9/10 mm interference screw
B-II-0815	K-wire Ø0.8 mm, length 150 mm, for 4/5/6 mm interference screw
B-II-0815B	K-wire Ø0.8 mm, length 150 mm, blunt ends, for 4/5/6 mm interference screw
B-II-1624	K-wire Ø1.6 mm, length 240 mm, for 7/8/9/10 mm interference screw
B-II-1624B	K-wire Ø1.6 mm, length 240 mm, blunt ends, for 7/8/9/10 mm interference screw
B-II-4000	Drill bit Ø4.0 mm, cannulated, for 4 mm interference screws, quick coupling
B-II-5000	Drill bit Ø5.0 mm, cannulated, for 5 mm interference screws, quick coupling
B-II-6000	Drill bit Ø6.0 mm, cannulated, for 6 mm interference screws, quick coupling







Distributor:

**BlueRock Medical**  
Possenhofener Str. 5  
D-82340 Feldafing

☎ +49 (0) 8157 998 1238  
🖨 +49 (0) 8157 998 1242  
✉ info@bluerockmedical.de

**[www.bluerockmedical.de](http://www.bluerockmedical.de)**

**Hotline: 0800 588 8005**



Manufacturer:

**bioretec**

**Bioretec Ltd.**

Hermiankatu 22,  
FI-33720 Tampere, Finland